

# A Powerful Tool for the Study of Airway Biology: Human Tracheobronchial Epithelial Model Grown at the Air-liquid Interface

## Application Note

CORNING

*Xu Xuemei, Xiaojun Wang, and Liu Jian*  
Corning Incorporated  
Corning Life Sciences Asia Technology Center  
Shanghai, China

### Introduction

Normal human tracheobronchial epithelium protects the airways by providing a barrier against injury from external insults. Intact cellular organization and function of the epithelial barrier provides a sterile and stable inner environment for living organisms. Ciliated and nonciliated columnar cells, secretory cells, and basal cells are arranged in a pseudostratified pattern to form complex epithelium structure and function. Two distinct compartments, apical and basolateral are developed by all cells creating tight junctions and contacting basement membrane, while some cells will reach the luminal surface<sup>1</sup>. When exposed to particles such as pollutants, bacteria, or virus, these particles will be cleared by the co-operation of ciliated cells and secretory cells<sup>2</sup>.


Human tracheobronchial epithelial cells grown at the air-liquid interface (ALI) culture have emerged as a powerful tool for the study of airway biology. Air-liquid interface culture for 3 to 6 weeks can result in the formation of well-differentiated, polarized human tracheobronchial epithelial cultures that resemble *in vivo* pseudostratified mucociliary epithelium. *In vitro* models using primary cells and cell lines are essential for understanding the function of the epithelium. Studies with cell lines are convenient, reproducible, and inexpensive compared to experiments involving primary cells. Calu-3, an epithelial cell line derived from human lung adenocarcinoma, has been widely used in studying bronchial epithelium, its closest resemblance to this epithelium *in vivo*<sup>3</sup>. ALI models of culturing epithelial cells has shown its application in response to stimulations including wounding, anion secretion, and viral infection<sup>4,5,6</sup>.

The Transwell® permeable support insert from Corning Life Sciences is a useful tool that mimics an *in vivo* environment, which is why it has been widely used in a variety of studies, including drug transport/absorption, cell migration and differentiation, and co-culture applications. In addition, several *in vitro* models have also been constructed to better explore the natural physiology of cellular mechanism.

### Reagents and Materials

- ▶ Calu-3 cell line (ATCC Cat. No. HTB-55)
- ▶ MEM medium (Corning Cat. No. 10-010-CV)
- ▶ Fetal bovine serum (FBS, Corning Cat. No. 35-076-CV)
- ▶ Transwell permeable support insert (Corning Cat. No. 3470, 0.4 µm PET membrane)
- ▶ 24-well plate (Corning Cat. No. 3524)
- ▶ 96-well assay plate (Corning Cat. No. 3915)
- ▶ Lucifer yellow (Sigma Cat. No. L0144)
- ▶ HBSS (Corning Cat. No. 21-020-CV)
- ▶ HEPES (Corning Cat. No. 25-060-Cl)
- ▶ Voltohmmeter (MilliporeSigma)
- ▶ SpectraMax® M4 plate reader (Molecular Devices)

## Workflow



## Procedure

### 1. Model Generation

Cells were seeded into Transwell® inserts on day 0. Prior to seeding, Transwell inserts were pre-incubated at 37°C with culture medium for at least 30 minutes. Calu-3 cells were seeded into the apical chamber of the Transwell inserts at densities of  $0.25 \times 10^5$ ,  $0.5 \times 10^5$ ,  $1 \times 10^5$ ,  $2 \times 10^5$ , and  $4 \times 10^5$  cells/well with 150  $\mu$ L of culture medium. 800  $\mu$ L of medium was added to the basal chamber. The plate was placed in a humidified 37°C/5% CO<sub>2</sub> incubator and cultured overnight. On day 1 for the A-L interface, the culture medium in apical chamber was removed after seeding for 24 hours, and cells were cultured with 500  $\mu$ L culture medium in the basal chamber. For L-L interface, the medium in both apical and basal chambers were kept. Culture medium was changed every 2 to 3 days, and cultured for 3 weeks.

### 2. TEER Value Monitor

The TEER value was measured every 2 to 3 days using a Voltohmmeter (MilliporeSigma) to confirm the development of tight junctions. Both the basolateral and apical chamber were filled with fresh pre-equilibrated medium and TEER was read after allowing the cell culture to get stable potential for approximately

15 to 30 min. to avoid any temperature fluctuation-induced TEER changes, because TEER measurement is sensitive to temperature<sup>7</sup>. Medium was removed from the apical side of the insert immediately after TEER readings to re-establish the ALI condition.

### 3. Lucifer Yellow Permeability Assay

On day 21, the Lucifer yellow permeability assay was performed. To prepare the donor solution, 100 mM stock solution of Lucifer yellow was diluted with HBSS/10 mM HEPES buffer to 300  $\mu$ M. The receiver solution used was HBSS/HEPES buffer with a DMSO content matched to that in the donor solution. The Transwell inserts were washed once with HBSS/HEPES buffer, transferred into a new 24-well receiver plate with receiver solution, and then 50  $\mu$ L of the donor solution was added to the apical chamber. The plate was placed at 37°C with orbital shaking for 60 min. Then 100  $\mu$ L of the solution in the receiver plate was transferred to a black assay plate and the fluorescence was read at 480 nm/530 nm using a SpectraMax® M4 plate reader. In the meantime, a series of Lucifer yellow solution with different concentrations were diluted and the standard curve was plotted. The concentration of Lucifer yellow in the receiver plate was calculated from the standard curve.

#### 4. Cell Staining

After culturing in the Transwell® inserts, Calu-3 cells, together with the Transwell membrane, were cut from the insert and fixed in 10% formalin overnight at 4°C. After dehydration and paraffin embedment, 5 µm sections were prepared and stained with hematoxylin and eosin (H&E), periodic acid–Schiff (PAS) stain for neutral proteoglycans, or Alcian blue for acid proteoglycans, and viewed with a lighted microscope.

To demonstrate mucus cell differentiation, confluent layers of cells grown at ALI model were subjected to PAS staining to visualize neutral and acidic glycoproteins (likely mucins).

#### 5. Calculations

##### Transepithelial Electrical Resistance (TEER):

TEER of insert and medium alone was subtracted from measured TEER and calculated multiplying by the insert area.

$$TEER = (R_m - R_b) \times A$$

Where:

TEER is the transepithelial electrical resistance

$R_m$  is the resistance reading obtained for the cell monolayer

$R_b$  is the resistance reading obtained for the blank insert (without cells)

A is the surface area of insert filter membrane [cm<sup>2</sup>]

##### Apparent Permeability ( $P_{app}$ ):

$$P_{app}[\text{cm/s}] = \frac{\text{Flux} \times V_d}{t \times A} = \frac{dQ}{dt} \times \frac{1}{A \cdot C_0}$$

Where:

Flux is the fraction of the donated amount recovered in the receiver chamber

$V_d$  is the volume in the donor chamber [cm<sup>3</sup>]

$C_0$  is the initial concentration in the donor solution [mM]

A is the surface area of insert filter membrane [cm<sup>2</sup>]

t is the incubation time [sec.]

$dQ/dt$  is the amount of drug transported within a given time period [pmol/sec.]


#### Results

##### 1. Microscopic Observation


In the L-L interface model, Calu-3 cells became confluent and lacked clear structure after 3 weeks of culturing in the Transwell inserts. In the A-L model, protuberances were observed after 3 weeks cultivation (Figure 1).

##### 2. Evaluation of Cell Monolayer Integrity

Transepithelial electrical resistance (TEER) is a widely accepted quantitative technique to measure the integrity of tight junction dynamics in cell culture models of epithelial monolayers. TEER



**Figure 1. Calu-3 bronchial epithelial model.** Cultivation of Calu-3 cell at L-L model (upper) and A-L model (below). Microscopy images photos were taken at day 0, 1 week, 2 weeks, and 3 weeks under microscope.



**Figure 2.** TEER measurements of Calu-3 in L-L and A-L models. Time-course of TEER development for the L-L model (left) and A-L model (right) at day 7, 14, and 21. Seeding densities of Calu-3 cells:  $4e^5$  (red squares),  $2e^5$  (blue dots), and  $1e^5$  (green triangles) per insert.

values are strong indicators of the integrity of the cellular barriers before they are evaluated for transport of drugs or chemicals. TEER measurements can be performed in real-time without cell damage and generally are based on measuring ohmic resistance or measuring impedance across a wide spectrum of frequencies.


In the L-L interface model, the TEER values of Calu-3 cells reached more than  $1,000 \Omega \cdot \text{cm}^2$  at day 14 and remained stable until day 21. In the A-L interface model, high TEER values were reached in one week, then decreased and stabilized at approximately  $300 \Omega \cdot \text{cm}^2$  (Figure 2).

Lucifer yellow was also used as a paracellular marker for monolayer integrity to confirm the results of TEER measurements. Lucifer yellow travels across the cell monolayer only through passive paracellular diffusion and has low permeability. It is not able to pass across cell monolayers when tight junctions between cells are maintained. Figure 3 shows the permeability ( $P_{app}$ ) of Lucifer yellow for Calu-3 epithelial L-L and A-L models.  $P_{app}$  was much lower than  $1 \times 10^{-6} \text{ cm/s}$  for both L-L and A-L models. Tight junctions were formed after 3 weeks of cultivation in the A-L model.

### 3. Haematoxylin-eosin (HE) Staining

The results of HE staining indicated that pseudostratified columnar epithelium was formed in 3 weeks in the A-L model (Figure 4).


Cells exhibited poor morphological quality in the L-L model. They appeared dedifferentiated, flat, or dwarf, with reduced numbers of cilia and appearance of secretory protrusions in the L-L model; While in the A-L model, cells were highly polarized, confluent, and retained a mix of ciliated and secretory cells.




**Figure 3.** Lucifer permeability of the Calu-3 epithelial models after 1 week and 3 weeks.

### 4. Periodic Acid Schiff (PAS) and Alcian Blue Staining of Proteoglycans/Mucus

Neutral proteoglycans were stained with PAS (red) and nuclei were stained with nuclear fast red (Figure 5, left panel). Acid proteoglycans were stained with Alcian blue (blue; Figure 5, right panel). The results showed that secretory products and vesicles were visible in cells cultured at the A-L interface for 3 weeks.



**Figure 4.** Haematoxylin-eosin stain of Calu-3 cells at the L-L and A-L interfaces after 1 week and 3 weeks in culture. Haematoxylin-eosin stain at the L-L interface revealed a flat cell monolayer (A, B). At the A-L interface, pseudostratified mucociliary epithelium morphology was formed (C, D).



**Figure 5.** PAS (left) and Alcian blue stain (right) for Calu-3 cells at the A-L interfaces at 3 weeks in culture.

## Conclusion

The study of airway biology can be achieved by cultivation of human tracheobronchial epithelial cells, primary cells, and cell lines at the air-liquid interface. In the present study, a successful bronchial epithelial model had been set up by using Calu-3 cells cultivated at the air-liquid interface for 3 weeks. Several methods were performed to verify the integrity of this model, including measurement of TEER value for tight junctions between cell monolayers and Lucifer yellow permeability. TEER reflects the ionic conductance of the paracellular pathway in the epithelial monolayer. This method is non-invasive and can be applied to monitor liver cells during their various stages of growth and differentiation. However, various factors have been shown to affect TEER values, including temperature, medium formulation, passage number of cells, and usage of electrodes<sup>7</sup>. Histochemical staining confirmed the pseudostratified mucociliary epithelium morphology.

## References

1. Pezzulo AA, Starner TD, Scheetz TE, Traver GL, Tilley AE, Harvey BG, Crystal RG, McCary PB, Zabner J, The air-liquid interface and the use of primary cell cultures are important to recapitulate the transcriptional profile of in vivo airway epithelia. *Am J Physiol Lung Cell Mol Physiol*, 2011, 300:L25-L31.
2. Houtmeyers E, Gosselink R, Gayan-Ramirez G, Decramer M. Regulation of mucociliary clearance in health and disease. *Eur Respir J*, 1999, 13(5):1177.
3. Kreft ME, Jerman UD, Lasic E, Hevir-Kene N, Rizner TL, Peternel L, Kristan K, The characterization of the human cell line Calu-3 under different culture conditions and its use as an optimized in vitro model to investigate bronchial epithelial function, 2015, *Eur J Pharm Sci*, 69:1-9.
4. Devor DC, Singh AK, Lamber LC, DeLuca A, Frizzell RA, Bridges RJ, Bicarbonate and chloride secretion in Calu-3 human airway epithelial cells, *J General Physiol*, 1999, 113:734-760.
5. Harcourt JL, Caidi H, Anderson LJ, Hynes LM, Evaluation of the Calu-3 cell line as a model of in vitro respiratory syncytial virus infection, *J Virol Methods*, 2011, 174:144-149.
6. Zhu Y, Chidekel A, Shaffer TH, Cultured human airway epithelial cells (Calu-3): A model of human respiratory function, structure, and inflammatory responses, *Critical care research and practice*, 2010, 2010:394578.
7. Srinivasan B, Kolli AR, Esch MB, Abaci HE, Shuler ML, Hickman JJ, 2015, *J Lab Autom*, 20(2):107-126.

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

# CORNING

**Corning Incorporated**  
*Life Sciences*

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

### ASIA/PACIFIC

**Australia/New Zealand**  
t 61 427286832

**China**  
t 86 21 3338 4338  
f 86 21 3338 4300

**India**  
t 91 124 4604000  
f 91 124 4604099

### Japan

t 81 3-3586 1996  
f 81 3-3586 1291

### Korea

t 82 2-796-9500  
f 82 2-796-9300

### Singapore

t 65 6572-9740  
f 65 6861-2913

### Taiwan

t 886 2-2716-0338  
f 886 2-2516-7500

### EUROPE

[CSEurope@corning.com](mailto:CSEurope@corning.com)

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 020 655 79 28  
f 020 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

### All Other European Countries

t +31 (0) 206 59 60 51  
f +31 (0) 206 59 76 73

### LATIN AMERICA

[grupoLA@corning.com](mailto:grupoLA@corning.com)

**Brasil**  
t 55 (11) 3089-7400

**Mexico**  
t (52-81) 8158-8400