

Corning® Lambda™ EliteMax Semi-automated Benchtop Pipettor

Frequently Asked Questions

The logo consists of the word "CORNING" in a white, uppercase, sans-serif font, centered within a solid orange square.

General

1. How do I operate this instrument?

For operating instructions, refer to the Corning Lambda EliteMax Semi-automated Benchtop Pipettor Instruction Manual (CLS-AN-671DOC) at www.corning.com/lifesciences.

2. What items come with this instrument?

Each instrument is delivered with a touch screen controller, external power supply, portrait adapter plate, instruction manual, and calibration certificate. Note that due to the length of the instruction manual, it will come on a thumb drive in an effort to reduce printing and increase sustainability.

3. What plug types come with the external power supply?

Each instrument is delivered with the following plug types: United States (US), United Kingdom (UK), Europe (EU), and China (CN) plug types.

4. What are the minimum and maximum volumes the liquid handler dispenses?

The dispenser has a volume range from 1 to 200 μL . The unit has been tested down to 2 μL with accuracy and precision less than 5%, and at 1 μL has about 10% accuracy and precision.

5. What types of protocols is the unit capable of running?

Plate fillings, serial dilutions, plate-to-plate transfers, plate replication, sample preparation, cherry picking protocol (for single-channel head only), and other more complex protocols can be implemented. For performance specifications, refer to the Corning Lambda EliteMax Semi-automated Benchtop Pipettor Instruction Manual (CLS-AN-671DOC) at www.corning.com/lifesciences.

6. What is the maximum volume of tips for a mixing protocol?

When doing multiple aspirate and dispense, we recommend using 150 μL as the maximum since there is residual volume after each aspirate/dispense cycle.

7. What is the lowest volume I could use for mixing?

The recommended minimum volume for mixing is 5 μL . Depending on the volume in the plate, a volume lower than 5 μL may result in insufficient mixing.

8. Is the tip removal manual or can it be done using the software?

Manual ejection is not required as the tips are ejected automatically back into the tip rack.

9. Is it possible to pipet different volumes from different wells in the same plate?

This is possible when using the cherry picking protocol (for single-channel head only), but is not possible in other protocols.

Software

1. How do I access the software for the instrument?

The touch screen controller has the software preloaded on the device and no download is required.

2. What happens if there is a software update?

To be notified of any software updates, follow the registration instructions on the instruction manual thumb drive. For software update instructions, refer to the Corning® Lambda™ EliteMax Semi-automated Benchtop Pipettor Instruction Manual (CLS-AN-671DOC) at www.corning.com/lifesciences.

3. What do I do if the system crashes?

The instrument has collision detection features and will notify the user about a possible crash. For more details about the crash recovery routine, refer to the Corning Lambda EliteMax Semi-automated Benchtop Pipettor Instruction Manual (CLS-AN-671DOC) at www.corning.com/lifesciences.

Pipet Tips

1. What kinds of pipet tips are compatible with this instrument?

This instrument is designed to be used with Axygen® FX tips only and should not be operated with any other types of tips.

2. What are the Axygen FX tip types and tip size ranges?

The Axygen FX tips are polypropylene, 96-well ANSI/SLAS formatted racked tips with volumes of 20, 30, 50, 165, and 250 µL. Sterile and non-sterile versions are available for all tip sizes, and wide bore tips are available in both 165 and 250 µL sizes.

3. What types of reagents are compatible with the Axygen FX tips?

For a list of reagents compatible with polypropylene tips, refer to the Chemical Resistance Chart (CLS-DL-CC-082) at www.corning.com/lifesciences.

4. Can the tips be discarded after use?

Yes, tips can be discarded into an empty tip rack in the Serial Dilution protocol and when using a single-channel head. The single-channel can discard tips into an empty rack, whereas while using the Serial Dilution protocol the tips can be discarded in a waste reservoir in position AA, BB, CC, DD or EE, depending on the program.

Other

1. Can this instrument fit under a hood?

Yes. The instrument is designed to fit under any standard sized hood. The dimensions of the instrument and touch screen controller can be found in the Corning Lambda EliteMax Semi-automated Benchtop Pipettor Instruction Manual (CLS-AN-671DOC) at www.corning.com/lifesciences.

2. What plates and reservoirs are compatible with this instrument?

All standard ANSI/SLAS formatted plates and reservoirs are compatible with this instrument.

3. How do I switch from portrait to landscape position on the deck plate?

Use the adapter plate to transition between portrait and landscape positions. For more information, refer to the Corning Lambda EliteMax Semi-automated Benchtop Pipettor Instruction Manual (CLS-AN-671DOC) at www.corning.com/lifesciences.

4. Can this instrument be used with UV light?

Yes. The instrument can be used with UV light. Note that long exposure to UV light might fade the colors but won't affect the operation or performance of the instrument.

5. What is ANSI/SLAS format?

This standard defines the well center positional requirements of a microplate as specified in American National Standards covering these microplates (as referenced at www.slas.org).

6. What are the performance specifications and technical data?

For performance specifications and technical data, refer to the Corning® Lambda™ EliteMax Semi-automated Benchtop Pipettor Instruction Manual (CLS-AN-671DOC) at www.corning.com/lifesciences.

7. What happens when the power shuts off during a run?

When the power goes off during a run, the run is stopped and the progress is lost. When the power comes on, the unit will return to home the axis.

8. How can I change the plates from landscape position to portrait position?

The plate orientation is predetermined before running the protocol – either landscape or portrait. The portrait adapter is included with the instrument and goes in position BB when running portrait protocols.

9. How long does it take to seed the cells per plate?

The speed depends on the protocol, however it takes about 1.5 minutes to fill the entire plate.

10. Is it possible to remove and store the media supernatant from the spheroid plate on another plate?

Yes, a plate-to-plate transfer can achieve this. The spheroid plate should be set as the source plate, and the empty plate will be the destination plate.

11. Is the Corning EliteMax benchtop pipettor suitable for biological (human) plasma and serum?

The instrument has been confirmed to be able to handle cell culture which contains serum and should be suitable for working with plasma and serum.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA

t 800.492.1110
t 978.442.2200

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

Chinese Mainland
t 86 21 3338 4338

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CSEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400